Medical Isolator
MI 1005

Medical Isolator
MI 2005

Handbuch/ User Manual

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MI 1005

MI 2005

© Baaske Medical GmbH & Co. KG. All rights reserved.
Information in this document has been carefully checked for accuracy; however, no guarantee is given to the correctness of the contents. Photographs and schematic diagrams used in this document are for illustration purposes only. The actual product may vary. This document contains proprietary information protected by copyright. No part of this manual may be reproduced by any mechanical, electronic, or other means, in any form, without prior written permission of Baaske Medical.

This document is subject to change without notice.

Content

1. General Safety Advice	23
1.1 Place of installation.....	23
1.2 Connection and Operation.....	24
1.3 Used symbols	25
1.4 Conformity.....	26
2. Intended Use	27
3. Installation	27
4. Maintenance	28
5. Cleaning	29
6. Potential Risk and Errors	29
7. Technical Specification	30
8. Scope of delivery.....	34
9. Warranty Regulations	34
10. Manufacturer Contact	35

1. General Safety Advice

Read this manual carefully before installing and implementing the MI 1005/MI 2005.

- a. The MI 1005/MI 2005 must be installed and implemented by qualified and trained personnel.
- b. Do not open the case of the MI 1005/MI 2005. There are no serviceable parts inside.
- c. Faulty or defective devices have to be exchanged and returned to the manufacturer.
- d. Before operating the MI 1005/MI 2005 in combination with medical devices, make sure that the manufacturer of the medical device permits the galvanic isolation of the Ethernet Signal.

1.1 Place of Installation

During installation ensure that the MI 1005/MI 2005 is protected from penetration by moisture, high humidity and mechanical damage. It should not be operated in areas or atmospheres with combustible gas mixtures.

1. General Safety Advice

The MI 1005/MI 2005 should also not be operated in close proximity to devices that emit electromagnetically radiation to avoid electromagnetically interferences like for example mobile phones. Potential interferences can be avoided by operating the MI 1005/MI 2005 outside this electromagnetically environment.

1.2 Connection and Operation

When connecting the MI 1005/MI 2005 Medical Isolator make sure that:

- a. Not any telecommunication appliances, voltage supplies or similar are connected to the MI 1005/MI 2005. By doing so, both the MI 1005/MI 2005 and personnel could be harmed!
- b. All connected Ethernet components comply with IEEE 802.3 or IEEE 802.3ab (10/100/1000-BaseT, Twisted-Pair) or 802.3bz 2.5GBase-T (MI 2005)
- c. Tested and approved CAT 6 network cable (with appropriate test mark) is used.

1.3 Used symbols

	CE Mark. The Product complies with all applicable EU Directives
	UL Mark „Recognized Component“ for US and Canadian Market
	Environmental advise. Do not dispose the product in the garbage. Return the product to the manufacturer for disposal.
	RoHS Mark. The Product complies with all applicable EU RoHS Directives.
	Please note the accompanying documents!

1. General Safety Advice

1.4 Conformity

Your device meets the electromagnetic compatibility requirements and has been tested to the standards specified herein.

UL **AAMI** ES60601-1:2005,ES60601-1:2005/
AMD1 1:2012 ,
ES60601-1:2005/AMD2:2021,CAN/CSA-
C22.2 No. 60601-1:08,
CAN/CSA-C22.2 No. 60601-1:14 (including
amendment 1)
and Amendment 2:2022 (MOD) to CAN/
CSA-C22.2 No. 60601-1:14

TÜV **IEC 60601-1:2005/AMD1:2012/AMD2:2020**
SÜD (Ed. 3.2)

EMV **EN 60601-1-2:2015** (IEC 60601-1-2:2014)
including tests according to:
(IEC 61000-4-2:2008: IEC 61000-4-3:2006 +
A1:2007 + A2:2010, IEC 61000-4-4:2012,
IEC 61000-4-6:2013)
EN 55032:2015

Full declarations of conformity are available on
request.

Hotline +49 (0) 5741 2360270

2. Intended Use

The MI 1005/MI 2005 is an appliance especially developed for use in medical technology. The purpose of the MI 1005/MI 2005 is to connect, in compliance with the relevant norms for electrical safety (IEC 60601-1) and EMC (IEC 60601-1-2), a medical appliance or system equipped with a network interface to a non-medical appliance or system (Single PC or local network). Examples for the use of the MI 1005/MI 2005:

- Connection to the practice or clinic networks with the PC of a prolonged EEG system
- Connection to an analysis computer in the doctor's or consulting room with an EMG measuring station

3. Installation

Only qualified and trained personnel should install the MI 1005/MI 2005. Please instruct your network systems administrator or your medicine technology department for this purpose. The installation of the MI 1005/MI 2005 is done by means of the two pluggable Ethernet Interfaces (RJ-45-sockets). To do this, plug the network connection cable into the socket of the MI 1005/MI 2005, until the lock

3. Installation

engages with an audible ‚click‘ (**observe installation direction for MI 2005!**). The installation of the MI 1005/MI 2005 should be done near the medical appliance/system. The safety advice in Chapter-1.2 – Connection and Operation- is to be observed!

4. Maintenance

For maintaining a safe technical condition, an annual inspection is recommended, which includes the following checks:

1. Examination for external damage (housing, network connections, legibility of labeling, dirt, etc.) of the MI 1005/MI 2005.
2. Examination for penetrated substances, especially liquids.
3. Function test.
4. Isolation test. The isolation resistance of the MI 1005/MI 2005 should not exceed 20 MOhm.
For testing, at both connections all data cables as well as the screen should always be shut down briefly. The isolation resistance should be measured between signal entry and signal exit.
5. Availability and sufficiency of the documents.

4. Maintenance

Should any deficiency be established during a test, the MI 1005/MI 2005 should not be used any longer, as the safety of the patients and the operator cannot be guaranteed.

In case of doubt please contact your supplier or the manufacturer.

5. Cleaning


Use a dry or a moderately moisturized cloth to clean the MI 1005/MI 2005. It must be noted that no fluids or moisture will penetrate into the device.

6. Potential Risk and Errors

Not any Risk or Errors are expected when the MI 1005/MI 2005 is being operated in compliance with the instructions in this manual. The fault / error rate of the MI 1005/MI 2005 is lower than 1:100000.


7. Technical Specifications MI 1005

Mark of origin (Manufacturer)		Baaske Medical GmbH & Co. KG
Type		MI 1005
Operating Mode		Continuous Operation
Environmental Conditions		
During Operation	Temperature	-10°C up to +85°C
	Rel. Humidity	10% - 90% (not condensing)
	Air pressure	700 hPa up to 1060 hPa
During Storage / Transport	Temperature	-10°C up to +85°C
	Rel. Humidity	10% - 95%(not condensing)
	Air pressure	500 hPa up to 1060 hPa
Protection grade according EN 60529		IPX0
Dimensions (LDH)		65 x 29 x 23 mm
Weight		~ 30 g
Operation Voltage Max.		300 V AC / 400 V DC
Isolation Voltage AC 50Hz 60s		5,0 kV
Reinforced Isolation		according to IEC 60601-1

Supported Network protocols and transmission speed	10 Mbit/s 10Base-T (IEEE802.3 Cl.14) 100 Mbit/s 100Base-Tx (IEEE802.3 Cl.25) 1000 Mbit/s 1000Base-T (IEEE802.3 Cl.25)
Interfaces	2 x RJ45 port
Material	UL94V-0 complying materials RoHs compliant
MTBF (Mean time to failure)	0,21 x 10 ⁸ (1 Error in 21 000 000 hours)
Product Classification	Passive ethernet isolation device, Isolation rate DI (300V AC / 400V DC) tested with 5kV AC
ESD Safety (electrostatic discharge)	-
Means of Protection	2 MOPP
Isolation Diagram	 The diagram illustrates the electrical isolation between two Ethernet ports. It features a transformer symbol (two vertical lines of unequal length) enclosed within a rectangular box labeled 'Plastic Enclosure'. The left side of the box is labeled 'RMS Ethernet Input / Output' and the right side is also labeled 'RMS Ethernet Input / Output'. Above the transformer, a horizontal double-headed arrow points left and right, with the letter 'A' positioned in the center of the arrow.

7. Technical Specifications MI 2005

Mark of origin (Manufacturer)		Baaske Medical GmbH & Co. KG
Type		MI 2005
Operating Mode		Continuous Operation
Environmental Conditions		
During Operation	Temperature	-10°C up to +85°C
	Rel. Humidity	10% - 90% (not condensing)
	Air pressure	700 hPa up to 1060 hPa
During Storage / Transport	Temperature	-10°C up to +85°C
	Rel. Humidity	10% - 95% (not condensing)
	Air pressure	500 hPa up to 1060 hPa
Protection grade according EN 60529		IPX0
Dimensions (LDH)		65 x 29 x 23 mm
Weight		~ 30 g
Operation Voltage Max.		300 V AC / 400 V DC
Isolation Voltage AC 50Hz 60s		5,0 kV
Reinforced Isolation		according to IEC 60601-1

Supported Network protocols and transmission speed	10 Mbit/s 10Base-T (IEEE802.3 Cl.14) 100 Mbit/s 100Base-Tx (IEEE802.3 Cl.25) 1000 Mbit/s 1000Base-T (IEEE802.3 Cl.25) 2500 Mbit/s 2.5GBase-T (IEEE802.3bz)
Interfaces	2 x RJ45 port
Material	UL94V-0 complying materials RoHs compliant
MTBF (Mean time to failure)	0,21 x 10 ⁸ (1 Error in 21 000 000 hours)
Product Classification	Passive ethernet isolation device, Isolation rate DI (300V AC / 400V DC) tested with 5kV AC
ESD Safety (electrostatic discharge)	Air 15 kV, Contact 8 kV
Means of Protection	2 MOPP
Isolation Diagram	 The diagram illustrates an isolation circuit for an Ethernet device. It features a transformer with two primary windings and one secondary winding. The primary windings are connected to two RJ45 Ethernet ports, labeled 'RJ45 Ethernet Output' and 'RJ45 Ethernet Input'. The secondary winding is connected to a central point, which is then connected to a 'Plastic Enclosure'. The circuit includes several diodes and a central component, likely a microcontroller or a relay, to manage the isolation and data flow between the two ports.

8. Scope of Delivery

- MI 1005/MI 2005 Network isolator
- User Manual
- connection cable CAT6 0,25m grey PIMF (optional)

Check the delivery contents and request any missing items from your retailer immediately.

9. Warranty Regulations

The manufacturer of the device is only responsible for damage in terms of safety, reliability and performance if repairs are carried out by persons authorized by him and the device is used in accordance with the instructions for use.

The terms and conditions of Baaske Medical GmbH & Co. KG, valid at the time of delivery of the product, are available under the Internet address:

www.baaske-medical.com/general-terms-and-conditions

10. Manufacturer Contact

For queries and questions please contact:

Manufacturer:

Baaske Medical GmbH & Co. KG

Bacmeisterstr. 3

32312 Lübbecke

Germany

Phone: +49 5741 236027- 0

Fax: +49 5741 236027- 99

E-mail: sales@baaske-medical.com

Web: www.baaske-medical.com

Baaske Medical GmbH & Co. KG

Bacmeisterstr. 3

32312 Lübbecke, Germany

Phone: +49 5741 236027- 0

Fax.: +49 5741 236027- 99

E-mail: sales@baaske-medical.com

www.baaske-medical.com